

# Freedom Excel Superlite Self-Propelled Wheelchair - 18"

SKU: MWS452170

**\$697**



## Features and Benefits:

- Folding backrest for easy transportation
- Flip up removable padded armrest
- Park brake below seat level when engaged for easy side transfer
- Outwards swinging removable legrest
- Height Adjustable Plastic Footplates
- Airless "soft ride" polyurethane tyres
- Quick release removable spoke wheels
- Double-braced aluminium construction
- Spill proof nylon fabric upholstery
- Footrest Length: 300-400mm
- Folded Width: 290 mm

## Specifications:

<b>Code</b>	MWS452170	<b>Safe Working Load</b>	110kg
<b>Description</b>	Freedom Excel Superlite Self-Propelled Wheelchair – 18" (46cm)	<b>Unit Weight</b>	12.7kg
<b>Total Length</b>	1070mm	<b>Seat Depth</b>	430mm
<b>Total Width</b>	640mm	<b>Seat Width</b>	460mm
<b>Total Height</b>	950mm	<b>Seat to Floor Height</b>	500mm
<b>Seat to Footrest</b>	500mm	<b>Backrest Height</b>	450mm