

Freedom Excel Superlite Transit Wheelchair - 20"

SKU: MWS452120

\$697



Features and Benefits:

- Folding backrest for easy transportation
- Handbrake with park locking mechanism connected to hub brakes
- Double-braced aluminium construction
- Drum hub brake system
- Height adjustable plastic footplates
- Outwards swinging removable legrest
- Airless "soft ride" polyurethane tyres
- Spill proof nylon fabric upholstery
- Footrest Length: 300-400 mm
- Folded Width: 290 mm

Specifications:

Code	MWS452110	Safe Working Load	110kg
Description	Freedom Excel Superlite Transit Wheelchair – 20" (51cm)	Unit Weight	11.8kg
Total Length	920mm	Seat Depth	430mm
Total Width	640mm	Seat Width	510mm
Total Height	950mm	Seat to Floor Height	500mm
Seat to Footrest	500mm	Backrest Height	450mm