

Freedom Premium Over Toilet Aid – HBA400

SKU: BTT163130

\$138



Features and Benefits:

- Splash guard included
- Slip resistant suction rubber tips
- Height adjustable legs
- Thicker, more durable powder coated aluminium frame to prevent scuffs and scratches
- Larger seat made from high density plastic to minimise flex; removable for easy cleaning

Specifications:

Code	BTT163130	Overall Height	650mm – 800mm
Description	Freedom Premium Over Toilet Aid	Seat Width	360mm
Overall Width	580mm	Seat Depth	410mm
Overall Depth	610mm	Seat Height	440mm – 590mm
Safe Working Load	140kg		