

# Freedom Premium Shower Chair – HBA402

SKU: BTS150110

**\$145**



## **Features and Benefits:**

- Solid plastic seat, perforated for drainage
- Thicker, more durable powder coated aluminium frame to prevent scuffs and scratches
- S-bend shaped frame, provides maximum space between the armrests
- Slip resistant suction rubber tips
- Height adjustable legs

## **Specifications:**

<b>Code</b>	BTS150110	<b>Seat Width</b>	490mm
<b>Description</b>	Freedom Premium Shower Chair	<b>Seat Depth</b>	410mm
<b>Total Length</b>	640mm	<b>Seat to Floor Height</b>	440mm – 590mm
<b>Total Width</b>	600mm	<b>Safe Working Load</b>	140kg
<b>Total Height</b>	690mm – 940mm	<b>Unit Weight</b>	3.3kg