

# Aspire Raphael Quattro Lift Chair Four Action - Large

SKU: CHP227310 (Bronze)

SKU: CHP227311 (Marina Grey)



**\$3,567**

## Features and Benefits:

- Four motor lift recline chair for infinite positioning.
- Microsuede, leather-look upholstery
- Backlit hand control with integrated USB charging point
- Independent adjustable headrest, back, leg & lumbar support functions (4 Motors)
- Integrated lift & recline functions for easy sit-to-stand & maximum comfort.
- Intuitive headrest support, ideal when watching TV when reclined.
- Head & armrest covers for improved longevity.
- Part of the Aspire GOLD range of lift recline chair.
- Available in Bronze (Microsuede Leather-look), or Marina Grey (Microsuede Leather-look)



Bronze  
Microsuede



Marina Grey  
Microsuede

## Specifications:

<b>Code</b>	CHP227310 (Bronze) CHP227311 (Marina Grey) CHP227312 (Nordic Stone)	<b>Seat Width</b> <b>Seat Depth</b>	550mm 530mm
<b>Description</b>	Aspire Raphael Quattro Lift Chair – Medium – Four Action	<b>Backrest Height</b>	700mm
<b>Total Width</b>	810mm	<b>Seat to Floor Height</b>	500mm
<b>Total Height</b>	1200mm	<b>Safe Working Load</b>	130kg
	450mm (Bronze) 480mm (Marina Grey)	<b>Colour</b>	Bronze Microsuede Marina Grey Microsuede