

Aspire Raphael Quattro Lift Chair Four Action - Large

SKU: CHP227310 (Bronze)

SKU: CHP227311 (Marina Grey)

SKU: CHP227312 (Nordic Stone)



\$3,270



Bronze
Microsuède



Marina Grey
Microsuède

Features and Benefits:

- Four motor lift recline chair for infinite positioning.
- Microsuède, leather-look upholstery
- Backlit hand control with integrated USB charging point
- Independent adjustable headrest, back, leg & lumbar support functions (4 Motors)
- Integrated lift & recline functions for easy sit-to-stand & maximum comfort.
- Intuitive headrest support, ideal when watching TV when reclined.
- Head & armrest covers for improved longevity.
- Part of the Aspire GOLD range of lift recline chair.
- Available in Bronze (Microsuède Leather-look), or Marina Grey (Microsuède Leather-look)

Specifications:

Code	CHP227310 (Bronze) CHP227311 (Marina Grey) CHP227312 (Nordic Stone)	Seat Width Seat Depth	550mm 530mm
Description	Aspire Raphael Quattro Lift Chair – Medium – Four Action	Backrest Height	700mm
Total Width	810mm	Seat to Floor Height	500mm
Total Height	1200mm	Safe Working Load	130kg

	450mm (Bronze) 480mm (Marina Grey)	Colour	Bronze Microsuede Marina Grey Microsuede
--	---------------------------------------	---------------	---