

Aspire Raphael Quattro Lift Chair Four Action - Medium

SKU: CHP227300 (Bronze)

SKU: CHP227301 (Marina Grey)

SKU: CHP227302 (Nordic Stone)



\$3,195

Features and Benefits:

- Four motor lift recline chair for infinite positioning.
- Microsuede, leather-look upholstery
- Backlit hand control with integrated USB charging point
- Independent adjustable headrest, back, leg & lumbar support functions (4 Motors)
- Integrated lift & recline functions for easy sit-to-stand & maximum comfort.
- Intuitive headrest support, ideal when watching TV when reclined.
- Head & armrest covers for improved longevity.
- Part of the Aspire GOLD range of lift recline chair.
- Available in Bronze (Microsuede Leather-look), or Marina Grey (Microsuede Leather-look)



Bronze
Microsuede



Marina Grey
Microsuede

Specifications:

| | | | |
|---------------------|---|--|----------------|
| Code | CHP227300 (Bronze) CHP227301 (Marina Grey) CHP227302 (Nordic Stone) | Seat Width Seat Depth | 520mm 500mm |
| Description | Aspire Raphael Quattro Lift Chair – Medium – Four Action | Backrest Height | 810mm |
| Total Width | 810mm | Seat to Floor Height | 470mm |
| Total Height | 1150mm | Safe Working Load | 130kg |

| | | | |
|--|---------------------------------------|---------------|---|
| | 450mm (Bronze) 480mm (Marina Grey) | Colour | Bronze Microsuede Marina Grey Microsuede |
|--|---------------------------------------|---------------|---|