

Aspire Raphael Quattro Lift Chair – Four Action - Small

SKU: CHP227290 (Bronze Microsuede)

SKU: CHP227291 (Marina Grey Microsuede)



\$3,487



Bronze Microsuede



Marina Grey Microsuede

Features and Benefits:

- Four motor lift recline chair for infinite positioning.
- Microsuede, leather-look upholstery
- Backlit hand control with integrated USB charging point
- Independent adjustable headrest, back, leg & lumbar support functions (4 Motors)
- Integrated lift & recline functions for easy sit-to-stand & maximum comfort.
- Intuitive headrest support, ideal when watching TV when reclined.
- Head & armrest covers for improved longevity.
- Part of the Aspire GOLD range of lift recline chair.
- Available in Bronze (Microsuede Leather-look), or Marina Grey (Microsuede Leather-look)

Specifications:

Code	CHP227290 (Bronze) CHP227291 (Marina Grey)	Seat Width	460mm (Bronze) 480mm (Grey)
Description	Aspire Raphael Quattro Lift Chair – Small – Four Action	Backrest Height	630mm
Total Width	740mm	Seat to Floor Height	470mm
Total Height	1030mm (Bronze) 1060mm (Marina Grey)	Safe Working Load	130kg
Seat Depth	450mm (Bronze) 480mm (Marina Grey)	Colour	Bronze Microsuede Marina Grey Microsuede