

# Aspire Raphael Quattro Lift Chair – Four Action - Small

SKU: CHP227290 (Bronze Microsuede)

SKU: CHP227291 (Marina Grey Microsuede)

## \$3,195



Bronze  
Microsuede



Marina Grey  
Microsuede

### Features and Benefits:

- Four motor lift recline chair for infinite positioning.
- Microsuede, leather-look upholstery
- Backlit hand control with integrated USB charging point
- Independent adjustable headrest, back, leg & lumbar support functions (4 Motors)
- Integrated lift & recline functions for easy sit-to-stand & maximum comfort.
- Intuitive headrest support, ideal when watching TV when reclined.
- Head & armrest covers for improved longevity.
- Part of the Aspire GOLD range of lift recline chair.
- Available in Bronze (Microsuede Leather-look), or Marina Grey (Microsuede Leather-look)

### Specifications:

<b>Code</b>	CHP227290 (Bronze) CHP227291 (Marina Grey)	<b>Seat Width</b>	460mm (Bronze) 480mm (Grey)
<b>Description</b>	Aspire Raphael Quattro Lift Chair – Small – Four Action	<b>Backrest Height</b>	630mm
<b>Total Width</b>	740mm	<b>Seat to Floor Height</b>	470mm
<b>Total Height</b>	1030mm (Bronze) 1060mm (Marina Grey)	<b>Safe Working Load</b>	130kg
<b>Seat Depth</b>	450mm (Bronze) 480mm (Marina Grey)	<b>Colour</b>	Bronze Microsuede Marina Grey Microsuede