

# Aspire Homecare Shower Chair

SKU: BTS118025

**\$115**



## **Features and Benefits:**

- Provides a complete seating support for users with reduced mobility during showering.
- Durable high grip feet for stability and safety.
- Adjustable legs to suit a wide range of users.
- Raised armrests assist the user in transferring on and off.
- Width between the arms: 465 mm.

## **Specifications:**

<b>Code</b>	BTS118025	<b>Seat Width</b>	410mm
<b>Description</b>	Aspire Homecare Shower Chair	<b>Seat Depth</b>	340mm
<b>Overall Width</b>	525mm	<b>Seat Height</b>	440mm – 570mm
<b>Overall Depth</b>	340mm	<b>Safe Working Load</b>	100kg
<b>Overall Height</b>	750mm – 880mm	<b>Width Between Arms</b>	465mm