

Pride Go-Go Boot Scooter

SKU: SCT670613

\$1,849



Features and Benefits

The Boot Scooter is loaded with design features like one-hand feather-touch disassembly and a convenient drop-in battery box that makes transport and travel worry free. The Boot Scooter delivers these features and more, along with the performance you expect from the first name in travel mobility, making it the ultimate travel scooter value.

- Feather touch disassembly into 5 pieces for easy transportation and storage (refer to image)
- Prides exclusive black non-scutting tyres
- Comfortable front-mounted post seat with arm rests for stability
- Easy to remove rechargeable battery
- Soft and easy to adjust tiller
- Flat free tyres
- Front basket standard
- Modular design for easy servicing
- Convenient off-board dual voltage charger that can charge battery pack on-board or off-board

Code	SCT670613	Safe Working Load	118kg
Description	Go-Go Boot Scooter – 4 Wheel	Maximum Speed	6.4km/hr
Total Length	1003mm	Range	11.5km
Total Width	495mm	Unit Weight	36.6kg (without batteries)
Colour	Red		