

Aspire Vogue Lightweight 2 Seat Walker

SKU: WAF705310CH (Champagne)

SKU: WAF705310RE (Red)

\$497



Features and Benefits

- Advanced locking breaks
- Unbreakable TPE super soft wheels
- New wider curved backrest
- Minimised storage space with folded width 220 mm
- Locks when folded
- Removable utility bag with 5 kg limit
- Cane holder attachment included

Specifications

Code	- WAF705310CH (Champagne) - WAF705310RE (Red)	Safe Working Load	150kg
Description	Vogue Lightweight 2	Unit Weight	6.1kg
Total Length	655mm	Seat to Floor Height	540mm
Total Width	620mm	Frame Weight	5.3kg
Total Height	785mm – 920mm	Seat Width	465mm

