

# Karma Ergo Lite Deluxe Self-Propelled Wheelchair – 18”

SKU: KCKM2512F20-18

**\$1,179**



**Features and Benefits:**

- Ultra lightweight, folding, attendant propelled or self-propelled style wheelchair
- S-Ergo patented seating system for added comfort
- Push and pull brake system for 14" rear wheel
- Detachable swing away footrests
- Optional fixed or quick release 20" self-propel wheel
- Armrests designed to allow closer seating to desks
- Top backrest section and frame folds for ease of transport
- S-Ergo patented seating system
- Adjusts and stabilises seating pose
- Double cross braces provide strength and durability

**Specifications:**

|                        |       |                             |        |
|------------------------|-------|-----------------------------|--------|
| <b>Max User Weight</b> | 100kg | <b>Seat Depth</b>           | 420mm  |
| <b>Overall Length</b>  | 930mm | <b>Floor to Seat Height</b> | 480mm  |
| <b>Overall Width</b>   | 645mm | <b>Folded Width</b>         | 295mm  |
| <b>Overall Height</b>  | 900mm | <b>Folded Depth</b>         | 695mm  |
| <b>Seat Width</b>      | 450mm | <b>Armrest Height</b>       | 215mm  |
| <b>Backrest Height</b> | 420mm | <b>Product Weight</b>       | 10.3kg |