

# Aspire Vogue Indoor Walker

SKU: WAF705100CH/RE

**\$425**



## Features and Benefits:

- **Luxury look and feel** – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- Now supplied with bonus cane holder (ref. image)
- **Lightweight and ergonomic design** with small turning circle for excellent manoeuvrability.
- **No lifting & spillage** – meals can be collected and brought to the desired location without the need for lifting or carrying. The tray is also removable.
- **Sitting down for a meal** – when the user is sitting, the walker can be turned around and positioned comfortably in front of them, to enjoy their meal.
- **Unique folding design** – the slimline design reduces the footprint when folded but allows the folded walker to stand upright on its own.
- **Brakes** – the ergonomic ‘Soft Touch’ lockable design allows the walker to be secured in place with minimal hand pressure.
- **Removable storage bag** – utensils and other items can be stowed securely away, ideal for transporting personal belongings around the house.

## Specifications:

<b>Total Length</b>	650mm	<b>Total Width</b>	540mm
<b>Total Height</b>	860mm – 970mm	<b>Colour</b>	Champagne / Red
<b>Safe Working Load</b>	110kg	<b>Unit Weight</b>	5.2kg