

# Aspire VIDA Folding Wheelchair – Self-Propelled

SKU: MWS449860 (Orange)

SKU: MWS449865 (Green)

**\$725**



## Features and Benefits:

- The Aspire Vida is a super lightweight transport wheelchair ideal for quick trips to the shops, appointments and catching up with friends.
- Ergonomic backrest and seat shape follow the contours of the user to maximise comfort and postural support.
- The 135kg safe working load caters to a broad range of users.
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers.
- Contoured and padded seating improves comfort and postural support.
- Attendant operated hand brakes with parking lock.
- Simple and intuitive folding for easy transportation.
- Removable leg rests and folding back for compact storage.
- Rear storage pocket for essentials.
- Available as attendant propelled or self-propelled.
- Anterior Seat Height - 495 mm
- Posterior Seat Height - 445 mm

## Specifications:

<b>Total Length</b>	1030mm	<b>Colour</b>	Orange, Green
<b>Total Width</b>	685mm	<b>Seat Depth</b>	420mm
<b>Total Height</b>	930mm	<b>Seat Width</b>	450mm
<b>SWL</b>	135kg	<b>Frame Weight</b>	10.7kg
<b>Unit Weight</b>	11.7kg	<b>Backrest Height</b>	430mm