

Aspire 46cm Weight Bearing Platform Shower Commode – Self-Propelled

Inc. Open Seat

SKU: BTC066200



\$2,175



Features and Benefits:

- 8° backrest angle ensures patient comfort
- Improved over toilet clearance
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- 2x Tente total locking castors and 2x large rear push wheels for self-propelling
- Safety arms prevent falls and swing back out of the way for transfers
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Easy-slide weight bearing platform footrest slides completely away for transfers
- Overall Height: 995 mm
- Overall Width: 665 mm
- Seat to Backrest: 305 mm
- Armrest Height: 225 mm
- Anterior Seat Height: 590 mm
- Posterior Seat Height: 560 mm

Specifications:

Code	BTC066200	Backrest Height	305mm
Overall Height	995mm	Armrest Height	225mm
Overall Width	665mm	Anterior Seat Height	590mm
Seat Width	460mm	Posterior Seat Height	560mm
Safe Working Load	200kg		