

Freedom Quad Adjustable Seat Walker

SKU: WAF750530 (Red)

\$237



Features and Benefits

- Seat walker / Rollator with adjustable seat height with 3 options 500, 540 or 580 mm
- Folding walker for easy transport and storage
- Ergonomic Handles
- Light touch lock brakes
- Foldable padded backrest
- Nylon storage bag
- Easy to fold without removing storage bag
- Height adjustable handles

Specifications:

Code	WAF750530 (red), BRO201-028	Seat to Floor Height	500mm – 580mm
Description	Freedom Quad Adjustable Seat Walker	Safe Working Load	140kg
Total Length	695mm – 780mm	Unit Weight	7.8kg
Total Width	620mm	Seat Width	440mm
Total Height	770mm – 960mm		