

# Days Whirl Transit Wheelchair – 18”

SKU: KUN-WHIRL45TR

**\$399**



## Features and Benefits:

- Has a range of features to accommodate different user needs.
- The Whirl features fully padded upholstery with a calf-support strap for comfort and safety.
- Features a detachable half-folding backrest to make reaching for objects behind you easier. The detachable backrest also makes the chair more compact for storage.
- Quick release wheels (for transport)
- Removable armrests
- Swing away and removable leg rests
- Solid tyres (no maintenance required)
- Front Castor Wheel: 8" x 1" PU
- Rear Wheel: 24" x 1 3/8" PU
- Upholstery Nylon (Flame Retardant)
- Colour: Silver
- Warranty: 12 months on frame & crossbars

## Specifications:

<b>Brand</b>	Days Healthcare	<b>Seat Width</b>	381mm
<b>Max SWL</b>	115kg	<b>Seat Depth</b>	430mm
<b>Seat Angle</b>	3°	<b>Backrest Height</b>	440mm
<b>Seat Height (Front)</b>	500mm	<b>Chair Weight</b>	14.5kg
<b>Seat Height (Rear)</b>	475mm		