

# Ellipse SuperLite Carbon Fibre Rollator - Small

SKU: 9346376000635 (Black)



**\$750**



## Features:

- ▶ Aust. STD Tested 492484
- ▶ Light Weight Carbon Fibre
- ▶ Added Cane Holder
- ▶ Easily Installable Nylon Bag
- ▶ Adjust Brake Tool-Free
- ▶ Adjust Height Tool-Free













## Accessories:

- ▶ Cup Holder- 9346376000574



- ▶ LED Light- LED-C-FX



Length (cm)	 66	Folded Width (cm)	 26	Front Castor Diameter (cm)	 20/8"
Width (cm)	 60	Seat Depth (cm)	 22	Rear Castor Diameter (cm)	 20/8"
Max Height (cm)	 85	Seat Width (cm)	 45	Max Load (kg)	 120
Minimum Height (cm)	 76	Seat Height from the Ground (cm)	 48	Net Weight (kg)	 6.26