

# Katie High Back Chair

SKU: RG635



**\$425**



- Features Include :**
- Modern design styling for a day chair
  - Adjustable seat height
  - Adjustable seat depth
  - Contoured padded backrest
  - Comfortable padded seat
  - Large padded armrests
  - Chair mobility wheels on rear legs as standard
  - Backrest folds down for transit + storage
  - Strong lightweight construction
  - Epoxy powder coated framework
  - Vinyl covered upholstery in Grey, Blue or Sage Green

Overall Height	Overall Width	Overall Length	Seat Depth	Seat Width	Seat Height	Backrest Height	SWL	Product Weight
1140mm	665mm	650mm	435mm - 475mm	480mm	450mm - 575mm	700mm	200kg	13kg

