

# Breezy BasiX 2 Lightweight Self-Propelled Wheelchair – 18”

SKU: 074103-005

**\$597**



## Features and Benefits:

Available in Self Propelled models, the Breezy BasiX2 is a lightweight wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look.

- Adapted seat heights accommodating for foot propulsion
- Armrests conveniently flip back, for easier transfers
- Easy accessible wheel locks prevent the wheelchair from rolling forward or backward
- Angle Adjustable Footplates
- Folding Frame for convenient transportation and storage
- Comes with free standard basic foam cushion
- Swing Away and Adjustable Anti-Tips
- Safe Working Load 125 kg
- Fixed Backrest

## Specifications:

<b>Max User Weight</b>	125kg	<b>Seat Width</b>	18”
<b>Seat Depth</b>	16”-18” (Adjustable)	<b>Seat Height</b>	19” (Adjustable 18”-20”)
<b>Back Height</b>	16”-18”	<b>Armrest Height</b>	8”-9”
<b>Unit Weight</b>	17.5kg		