

Breezy BasiX 2 Lightweight Self-Propelled Wheelchair – 20”

SKU: 074103-006

\$597



Features and Benefits:

Available in Self Propelled models, the Breezy BasiX2 is a lightweight wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look.

- Adapted seat heights accommodating for foot propulsion
- Armrests conveniently flip back, for easier transfers
- Easy accessible wheel locks prevent the wheelchair from rolling forward or backward
- Angle Adjustable Footplates
- Folding Frame for convenient transportation and storage
- Comes with free standard basic foam cushion
- Swing Away and Adjustable Anti-Tips
- Safe Working Load 125 kg
- Fixed Backrest

Specifications:

Max User Weight	125kg	Seat Width	20”
Seat Depth	16”-18” (Adjustable)	Seat Height	19” (Adjustable 18”-20”)
Back Height	16”-18”	Armrest Height	8”-9”
Unit Weight	17.5kg		