

# Freedom Toilet Pan Surround – HBA418

SKU: BTT165110

**\$80**



## Features and Benefits:

- Fits around the toilet pan
- Padded hand grip with curves for safety and comfort
- Height adjustable legs
- Top of frame is slightly angled to provide leverage when rising from the toilet
- Sturdy aluminium construction
- Slip resistant suction rubber tips

## Specifications:

<b>Code</b>	BTT165110	<b>Overall Height</b>	660mm – 785mm
<b>Description</b>	Freedom Toilet Pan Surround	<b>Unit Weight</b>	2kg
<b>Overall Width</b>	540mm	<b>Warranty</b>	1 Year
<b>Overall Depth</b>	470mm		