

REDGUM Space Saver Shower Chair (Caddy Included)

SKU: RG555H

\$149



Features and Benefits:

- Easy to assemble – NO tools required
- Adjustable in height with Sure Lock Adjustment
- Rubber leg ends for stability and grip
- Angled seat for maximum comfort
- Reinforced plastic seat
- Ribbed arm handles
- Strong lightweight construction
- Easy Clean Design
- Lifetime warranty against corrosion
- Seat and handles: Plastic
- Legs: Anodized Aluminium
- Leg Tips: Flexible PVC

Specifications:

Overall Dimensions (at lowest setting)	515mm (W) x 450mm (D) x 825mm (H)
Seat Depth	355mm
Seat Height	485mm – 610mm (from the floor to the top of the seat). Adjustable height in 25mm adjustment increments
Seat Handle	155mm (From top of seat to top of handle)
Product Weight	3.1kg
Maximum User Weight	130kg

