

# Pharmaquip Standard Lightweight Aluminium Self- Propelled Wheelchair

SKU: PA150

**\$415**



## Features and Benefits:

- Powder coated Aluminium wheelchair
- Double cross bar
- Plastic handgrips with attendant handbrakes
- Fixed armrest - PU arm pad
- Swing away detachable footrest
- Fold-down backrest
- Nylon upholstery
- PVC front castors
- Rear spoke wheel with solid PU tyre
- Conditional warranty applies

## Specifications:

<b>Seat Width</b>	46cm	<b>Folded Width</b>	29cm
<b>Seat Height</b>	48cm	<b>Overall Length</b>	104cm
<b>Seat Depth</b>	40cm	<b>Overall Height</b>	87.5cm (back folded: 70cm)
<b>Backrest Height</b>	40cm	<b>Front Castor</b>	20cm x 2.5cm
<b>Armrest Height</b>	20cm	<b>Rear Wheel Size</b>	60cm x 3cm
<b>Overall Width</b>	64.5cm	<b>Max Weight Capacity</b>	110kg
<b>Net Weight</b>	13.8kg	<b>Colours Available</b>	Silver or White Powder Coated