

# Top Gun Outlaw Scooter

SKU : OUTLAW

**\$3,470**



## Features

- Australia's Best Value!
- SWL – 150kg
- Mag. Wheels
- Front Lockbox
- Stylish!
- Available In Red or Blue

<b>Code</b>	OUTLAW	<b>Safe Working Load</b>	150kg
<b>Description</b>	Outlaw Scooter – 4 Wheel	<b>Maximum Speed</b>	10km/hr
<b>Total Length</b>	134.1cm	<b>Range</b>	35kms
<b>Total Width</b>	65cm	<b>Unit Weight</b>	106kg
<b>Colour</b>	Red or Blue	<b>Operation Time (when fully charged)</b>	Approx. 3.5 hours

