

# Folding Shower Chair

SKU: RG9400

**\$179**



## Features and Benefits:

The REDGUM Brand Foldable Shower Chair has been designed to easily fold for storage when not required.

This foldable shower chair has tough aluminium painted frame with stainless steel hinge pins with a polyurethane self draining seat.

The legs are easily adjusted to suit multiple seat height requirements and can be adjusted to accommodate uneven floors.

Each chair leg is fitted with durable rubber stoppers to minimise any movement.

The foldable shower chair is finished a pearlescent white finish.

## Specifications:

Weight Capacity	Seat Width	Adjustment Range	Colour
100kg	470mm	460mm – 530mm	Pearl White